Sistema de implantes CONELOG®

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Referencias

- (1) Semper W, Heberer S, Mehrhof J, Schink T, Nelson K. Effects of repeated manual disassembly and reassembly on the positional stability of various implant-abutment complexes: an experimental study. *Int J Oral Maxillofac Implants.* 2010;25(1):86-94.
- (2) Semper-Hogg W, Kraft S, Stiller S, Mehrhof J, Nelson K. Analytical and experimental position stability of the abutment in different dental implant systems with a conical implant-abutment connection. *Clin Oral Investig.* 2013;17(3):1017-23.
- (3) Semper-Hogg W, Zulauf K, Mehrhof J, Nelson K. The Influence of Torque Tightening on the Position Stability of the Abutment in Conical Implant-Abutment Connections. *Int J Prosthodont.* 2015;28(5):538-41.
- (4) Nelson K, Zabler S, Wiest W, Schmelzeisen R, Semper-Hogg W. Die Implantat-Abutment Verbindung. *Implantologie.* 2013;21(4):355-63.
- (5) Harder S, Dimaczek B, Acil Y, Terheyden H, Freitag-Wolf S, Kern M. Molecular leakage at implant-abutment connection--in vitro investigation of tightness of internal conical implant-abutment connections against endotoxin penetration. *Clin Oral Investig.* 2010;14(4):427-32.
- (6) Harder S, Quabius ES, Ossenkopf L, Kern M. Assessment of lipopolysaccharide microleakage at conical implant-abutment connections. *Clin Oral Investig.* 2012;16(5):1377-84.
- (7) Angermair J, Wiest W, Rack A, Zabler S, Fretwurst T, Nelson K. Synchrotron-based Radiography of Conical- vs. Butt-joint Implant Abutment Connections. *J Oral Implantol.* 2020 (online ahead of print).
- (8) Rack T, Zabler S, Rack A, Riesemeier H, Nelson K. An in vitro pilot study of abutment stability during loading in new and fatigue-loaded conical dental implants using synchrotron-based radiography. *Int J Oral Maxillofac Implants.* 2013;28(1):44-50.
- (9) Ackermann K-L, Barth T, Cacaci C, Kistler S, Schlee M, Stiller M. Clinical and patient-reported outcome of implant restorations with internal conical connection in daily dental practices: prospective observational multicenter trial with up to 7-year follow-up. *Int J Implant Dent.* 2020;6(1):14.
- (10) Lorenz J, Blume M, Korzinskas T, Ghanaati S, Sader RA. Short implants in the posterior maxilla to avoid sinus augmentation procedure: 5-year results from a retrospective cohort study. *Int J Implant Dent.* 2019;5(1):3.
- (11) Moergel M, Rocha S, Messias A, Nicolau P, Guerra F, Wagner W. Radiographic evaluation of conical tapered platform-switched implants in the posterior mandible: 1-year results of a two-center prospective study. *Clin Oral Implants Res.* 2016;27(6):686-93.
- (12) Moergel M, Nicolau P, Rocha S, Messias A, Guerra F, Wagner W. Clinical outcomes with CONELOG SCREW-LINE implants in the posterior mandible -final 5-year results of a prospective two-center study. *Int Poster J Dent Oral Med.* 2018;20(ORGs):Supplement, Poster 1213.
- (13) Messias A, Rocha S, Wagner W, Wiltfang J, Moergel M, Behrens E, et al. Peri-implant marginal bone loss reduction with platform-switching components: 5-Year post-loading results of an equivalence randomized clinical trial. *J Clin Periodontol.* 2019;46(6):678-87.
- (14) Schwarz F, Alcoforado G, Nelson K, Schaer A, Taylor T, Beuer F, et al. Impact of implant-abutment connection, positioning of the machined collar/microgap, and platform switching on crestal bone level changes. *Camlog Foundation Consensus Report.* *Clin Oral Implants Res.* 2014;25(11):1301-3.
- (15) Camlog. CAMLOG and Science. 2013;XJ6164.11/2013:22-32.
- (16) Cacaci C, Ackermann KL, Barth T, Kistler S, Stiller M, Schlee M. A non-interventional multicenter study to document the implants success and survival rates in daily dental practices of the CONELOG screw-line implant. *Clin Oral Investig.* 2019;23(6):2609-16.
- (17) Molina A, Sanz-Sánchez I, Martín C, Blanco J, Sanz M. The effect of one-time abutment placement on interproximal bone levels and peri-implant soft tissues: a prospective randomized clinical trial. *Clin Oral Implants Res.* 2017;28(4):443-52.
- (18) Fierravanti L, Ambrosio N, Molina A, Sanz I, Martín C, Blanco J, et al. The effect of one-time abutment placement on marginal bone levels and peri-implant soft tissues: 3 years results from a prospective randomized clinical trial. *Int Poster J Dent Oral Med.* 2018;20(ORGs):Supplement, Poster 1211.
- (19) Al-Sawaf O, Tuna T, Rittich A, Kern T, Wolfart S. Randomized clinical trial evaluating the effect of splinting crowns on short implants in the mandible 3 years after loading. *Clin Oral Implants Res.* 2020;00:1-11.
- (20) Schnutenhaus S, Brunkens L, Edelmann C, Dreyhaupt J, Rudolph H, Luthardt RG. Alveolar ridge preservation and primary stability as influencing factors on the transfer accuracy of static guided implant placement: a prospective clinical trial. *BMC Oral Health.* 2020;20(1):178.
- (21) Schnutenhaus S, Edelmann C, Rudolph H, Dreyhaupt J, Luthardt RG. 3D accuracy of implant positions in template-guided implant placement as a function of the remaining teeth and the surgical procedure: a retrospective study. *Clin Oral Investig.* 2018;22(6):2363-72.
- (22) Schnutenhaus S, Gröller S, Luthardt RG, Rudolph H. Accuracy of the match between cone beam computed tomography and model scan data in template-guided implant planning: A prospective controlled clinical study. *Clin Implant Dent Relat Res.* 2018;20(4):541-9.
- (23) Conserva E. Initial stability after placement of a new buttress-threaded implant. A case series study. *implants.* 2019(3):24-8.

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CAMLOG Biotechnologies GmbH | Margarethenstr. 38 | 4053 Basilea | Suiza
Fabricante de los productos CAMLOG® y CONELOG®: ALTATEC GmbH | Maybachstr. 5 | 71299 Winnenden | Alemania